



DISINTEGRATION TESTER

Disintegration Testers for Pharmaceutical Applications

A disintegration tester is a laboratory instrument used to determine how quickly tablets or capsules break down in a liquid medium. This is a critical quality control step for pharmaceutical manufacturers to ensure consistent drug release and proper absorption in the body.

Torontech's disintegration instrument offers reliable, precise testing with user-friendly controls and robust construction. It is suitable for routine testing in R&D, production, and QC labs.

Materials It Can Test

- Pharmaceutical tablets
- Hard and soft gelatin capsules
- Enteric-coated dosage forms
- Nutritional supplements
- Veterinary medications

Why Choose Torontech Disintegration Equipment

- **Accurate and Repeatable Testing:** Follows USP, EP, and other pharmacopoeial standards.
- **Independent Stations:** Multiple baskets allow parallel testing of different samples.
- **Automatic Time and Temperature Control:** Ensures stable conditions throughout the test cycle.
- **Easy Operation:** Digital display and simple controls for quick setup and monitoring.
- **Built-in Safety Features:** Overload protection and auto-stop improve lab safety.
- **Compact, Durable Design:** Space-saving construction ideal for busy laboratories.

Reliable Results for Consistent Quality

Torontech's disintegration tester helps pharmaceutical companies maintain strict quality control and regulatory compliance. With every test, the disintegration instrument provides fast, dependable insights into product performance.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

