

GUIDEWIRES TIP FLEXIBILITY TESTER



Guidewires Tip Flexibility Tester | ISO 11070, FDA Guidance (16007)

Guidewires Tip Flexibility Tester – ToronMed™ GTFT01 is a specialized instrument designed to evaluate the tip flexibility of intravascular guidewires and catheters. With its advanced displacement measurement system, high-precision load cell, and real-time data processing, this tester provides accurate and repeatable results for R&D laboratories, medical device manufacturers, and regulatory compliance testing.

Why Choose Guidewires Tip Flexibility Tester - ToronMed™ GTFT01?

The Guidewires Tip Flexibility Tester – ToronMed™ GTFT01 ensures precision, compliance, and reliability for evaluating guidewire and catheter tip flexibility. With its large touchscreen interface, automated calibration, transparent observation base, and GMP-ready data management, it is the ideal solution for medical device manufacturers, inspection authorities, and R&D centers. By ensuring flexibility standards are met, this tester helps guarantee safety, durability, and clinical performance of interventional devices.

APPLICATIONS

- Basic Applications
 - Intravascular guidewires
- Extended Applications
 - Thin catheters
 - Other interventional devices requiring tip flexibility testing



Standards Compliance

The ToronMed™ GTFT01 Guidewires Tip Flexibility Tester is designed in accordance with:

- ISO 11070 – Intravascular catheters – Sterile and single-use intravascular catheters
- FDA Guidance (16007) – Coronary, Peripheral, and Neurovascular Guidewires – Performance Tests and Recommended Labeling

FEATURES

- High Accuracy & Repeatability
Imported load cell ensures 0.5% FS accuracy, with stable performance and long service life.
- Advanced Displacement Control
Accuracy up to 0.1 mm, with stepless adjustment and precise control over the testing process.
- 10.1-Inch Color Touchscreen
Intuitive interface with real-time curve display and data analysis.
- Transparent Test Base
Professionally designed base allows clear observation of guidewire behavior during testing.
- Data Management
Stores 100,000+ test records, with micro printer and USB interface for convenient data transfer and reporting.
- User-Defined Testing
Adjustable curve points for customized test profiles and flexible research applications.
- Automated Calibration & Judgment
Simplifies operation while ensuring reliable performance.

- Data Security (Optional)
GMP-compliant software with user management, authority levels, and audit trail capabilities.
- Safety Protection
Over-range shutdown automatically prevents sensor overload.

THEORY & METHOD

Test Principle

The Guidewires Tip Flexibility Tester applies controlled force to the tip of a guidewire or catheter, measuring its displacement and resistance. This simulates clinical conditions to evaluate flexibility, durability, and overall product safety. Results are recorded, analyzed, and stored automatically.

TECHNICAL SPECIFICATIONS

Specification	Value
Test Range	0-500 g
Accuracy	0.5% FS
Resolution	0.001 g
Displacement Range	0.001-40.000 mm
Displacement Resolution	0.001 mm
Test Speed	0-10 mm/s

Specification	Value
Speed Accuracy	± 0.01 mm/s
Dimensions	500 × 400 × 350 mm
Net Weight	21 kg
Voltage	AC 220V 50Hz (110V is also available)

Standard Configuration

- Main instrument unit
- 10.1-inch touchscreen interface
- Transparent test base (2 pcs)
- Measuring block
- Micro printer
- Test fixture



TORONTECH

Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com



USA | Canada | EU | UAE | GCC | KSA | India | Africa | Latin America | Asia

