

BLOOD BAG EMPTYING UNDER PRESSURE TESTER



Blood Bag Emptying under Pressure Tester – Model ToronMed™ BBEUPT01 is designed to evaluate the pressurized emptying performance of plastic containers for single use, such as blood bags and related medical storage containers. By simulating real-world conditions, this tester ensures that blood bags and similar products maintain performance and safety during pressurized emptying, as required by ISO 3826-1 and other international standards.

Engineered with precision and reliability, the Blood Bag Emptying under Pressure Tester – Model ToronMed™ BBEUPT01 is an essential instrument for medical device manufacturers, laboratories, and regulatory facilities engaged in product validation, quality control, and compliance testing.

Torontech is a trusted global provider of advanced medical device testing solutions. With the Blood Bag Emptying under Pressure Tester – Model ToronMed™ BBEUPT01, laboratories and manufacturers gain a precise, durable, and standards-compliant solution that ensures the safety, performance, and reliability of single-use medical containers.

APPLICATIONS

The Blood Bag Emptying under Pressure Tester – Model ToronMed™ BBEUPT01 is ideal for:

- Medical device manufacturers producing blood bags and other single-use medical containers.
- Independent testing labs conducting compliance verification in line with ISO and international standards.
- Quality assurance departments performing routine validation and production line QC.
- Research institutions developing new blood bag materials or container designs.



FEATURES

- Pressurized Emptying Test – Evaluates container performance under controlled compression force and pressure.
- Adjustable Compression Speed – Provides a wide testing range to simulate various usage conditions.
- Accurate Pressure Measurement – High-resolution sensors ensure precise test data.
- Stable Environmental Performance – Operates reliably under controlled working temperatures.

TECHNICAL SPECIFICATIONS

Specification	Blood Bag Emptying under Pressure Tester - Model ToronMed™ BBEUPT01
Test Size	450 × 270 mm
Working Stroke	200 mm
Compression Speed	1-1000 mm/min
Pressure Loading Range	100 kPa
Pressure Resolution	0.1 kPa
Working Temperature	15-55 °C
Instrument Dimensions	Compression device: 500 × 500 × 750 mm Main control box: 300 × 380 × 150 mm
Net Weight	Compression device: 30 kg Main control box: 10 kg
Voltage	220V, 50Hz (110V is also available)
Standards	ISO 3826-1 and other applicable international standards.

Specification	Blood Bag Emptying under Pressure Tester - Model ToronMed™ BBEUPT01
Standard Configuration	Instrument, compression device.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

