



## DIGITAL WATER BATH

# Digital Water Bath – ToronWB™ D

## High-Precision Digital Water Bath for Demanding Laboratory Temperature Control

ToronWB™ D is a high-precision electric heating water bath instrument engineered for laboratory applications where tighter temperature control is required than standard water bath equipment can provide. With a temperature fluctuation of just  $\pm 0.2$  °C across its full RT+5 °C to 100 °C working range, the ToronWB™ D delivers the stability and reproducibility demanded by sensitive analytical, biological, and pharmaceutical procedures.

Constructed with a powder-coated steel housing and a corrosion-resistant stainless steel inner chamber, this water bath instrument combines robust build quality with the precise temperature regulation needed for exacting scientific work. Available in two sizes — the ToronWB™ D 500 and ToronWB™ D 700 — both models offer an open-top rectangular bath configuration that accommodates a wide variety of laboratory vessels, making this water bath an exceptionally versatile piece of equipment for any general-purpose or specialized laboratory.

## APPLICATIONS

---

### Digital Water Bath – ToronWB™ D Applications

- Precision enzyme kinetics and biochemical reaction incubation requiring tight temperature control
- Microbiological culture and pathogen growth studies at specific set-point temperatures
- Pharmaceutical dissolution testing and stability studies requiring  $\pm 0.2$  °C accuracy
- Molecular biology applications including hybridization, restriction digestion, and primer annealing
- Clinical laboratory serology and immunoassay plate incubation
- Chemical synthesis and reflux reactions requiring stable bath temperature
- Calibration support for thermometers and temperature-sensitive instruments

- Histology and cytology slide preparation and staining procedures



## FEATURES

---

### Digital Water Bath - ToronWB™ D Key Features

- Powder-coated quality steel housing provides corrosion resistance and durability in wet laboratory environments; stainless steel inner chamber ensures chemical resistance and ease of decontamination
- High-precision temperature sensor and integrated circuit control system deliver  $\pm 0.2$  °C temperature fluctuation — five times tighter than standard water bath equipment — for demanding analytical and biological applications
- Full operating range of RT+5 °C to 100 °C with 1 °C digital display resolution covers the complete spectrum of common laboratory heating requirements
- Open-top rectangular bath design accommodates a wide variety of laboratory vessels including test tubes, Erlenmeyer flasks, beakers, culture bottles, and multi-well plates
- Two model sizes (ToronWB™ D 500 and ToronWB™ D 700) provide flexibility to match bench space and sample throughput requirements
- Compact, easy-to-clean design with minimal internal cavities reduces contamination risk and simplifies routine maintenance

## THEORY & METHOD

---

### Theory and Method

ToronWB™ D uses indirect liquid-medium heating to deliver highly uniform, stable temperatures to submerged laboratory vessels. Distilled or deionized water fills the stainless steel inner chamber, which is heated by an immersion-type electric element. A high-precision RTD (Resistance Temperature Detector) or thermistor sensor monitors bath temperature continuously and relays real-time data to the integrated circuit control system.

Unlike standard water bath equipment with  $\pm 1\text{ }^{\circ}\text{C}$  tolerance, the ToronWB™ D series employs an enhanced control algorithm that modulates heater output with greater sensitivity, achieving temperature fluctuation of  $\pm 0.2\text{ }^{\circ}\text{C}$  — five times tighter than the base ToronWB™ EH series. This precision makes the ToronWB™ D water bath instrument the preferred choice for temperature-sensitive assays where small deviations can significantly affect experimental outcomes.

Temperature is set and displayed digitally with  $1\text{ }^{\circ}\text{C}$  resolution on the front-panel interface. The rectangular open-top bath design allows flexible positioning of tubes, flasks, beakers, and multi-well plates within the water bath, while the stainless steel chamber walls ensure even heat distribution throughout the bath volume. The quality steel housing with electrostatic powder-coat finish protects internal components from the humidity and chemical splashes inherent to typical water bath laboratory environments.

TECHNICAL SPECIFICATIONS

Digital Water Bath - ToronWB™ D Technical Specifications

| Parameter                   | ToronWB™ D 500                        | ToronWB™ D 700 |
|-----------------------------|---------------------------------------|----------------|
| Temp. Range                 | RT+5 - 100 °C                         |                |
| Temp. Resolution            | 1 °C                                  |                |
| Temp. Fluctuation           | $\pm 0.2\text{ }^{\circ}\text{C}$     |                |
| Power                       | 500 W                                 | 700 W          |
| Voltage                     | 220 V, 50 Hz (110V is also available) |                |
| Internal Dimensions (WxDxH) | 420x190x150 mm                        | 600x300x180 mm |
| Housing Dimensions (WxDxH)  | 565x220x260 mm                        | 745x340x295 mm |
| Net Weight                  | 8 kg                                  | 15 kg          |

| Parameter    | ToronWB™ D 500 | ToronWB™ D 700 |
|--------------|----------------|----------------|
| Gross Weight | 9 kg           | 16 kg          |
| Package Size | 670×330×340 mm | 860×450×390 mm |



**Torontech Inc.**

251 Consumers Road, Suite 1200  
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652  
Email: [info@torontech.com](mailto:info@torontech.com)

**Torontech, LLC**

601 Brickell Key Drive, Suite 700  
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919  
Email: [sales@torontech.com](mailto:sales@torontech.com)

**Torontech ME FZE**

Meydan Grandstand, 6th Floor,  
Meydan Road, Nad Al Sheba,  
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252  
Email: [me@torontech.com](mailto:me@torontech.com)

