



ANESTHESIA CART

Anesthesia Cart - ToronCAB™ 1300

ToronCAB™ 1300 hospital anesthesia cart organizes anesthesia drugs, airway supplies, and monitoring power for the operating room. Five color-coded drawers divide into sixteen adjustable compartments, so the team sorts drug classes by color and resizes each space to fit. An upper medicine shelf keeps ampoules, syringes, and supplies in clear bins above the work surface, while a side-pull work table adds preparation space. A power strip runs monitors at the cart, and four braking 125 mm casters move it and lock it steady.

Anesthesiologists and anesthesia technicians use the Torontech ToronCAB™ 1300 to keep drugs and equipment organized and within reach during a case. An impact-resistant HDPE body resists knocks and wipes down easily, while a welded steel internal frame keeps the cart rigid under heavy daily use. Drawers run on three-step, self-locking silent slides, and the side-pull table extends working space without enlarging the cart footprint. High-grade 125 mm silent casters let staff reposition the cart quietly in the operating room.

For North American facilities, a central control lock secures all five drawers in one action, supporting controlled-access handling of anesthesia drugs. Manufactured under ISO 13485 quality systems, the cart suits departments that want a certified, well-built anesthesia platform. The L815 by W540 mm footprint moves through standard doorways and positions easily within an operating room or procedure suite. Designed for daily clinical use, the ToronCAB™ 1300 brings classified drug storage, supply organization, power, and sharps disposal together on one station.

APPLICATIONS

Anesthesia Cart - ToronCAB™ 1300 Applications

The ToronCAB™ 1300 supports anesthesia and sedation workflows across departments where drugs, airway supplies, and monitoring all need to stay organized and mobile. Its color-coded drawers, medicine shelf, and power strip fit a range of

perioperative settings.

- Operating Rooms: keeps anesthesia drugs, syringes, and airway supplies organized and powered beside the anesthesia machine throughout surgery.
- Anesthesia Workrooms and Prep Areas: stocks and restocks the cart between cases, with color-coded drawers that speed setup for the next patient.
- Post-Anesthesia Care Units (PACU): brings drugs, monitoring power, and supplies to the bedside during recovery and emergence.
- Procedure and Sedation Suites: supports endoscopy, interventional, and imaging procedures that use sedation away from the main operating room.
- Intensive Care Units: provides organized drug storage and device power for sedation and airway management in critical care.
- Labor and Delivery: keeps anesthesia and analgesia supplies ready for epidurals, cesarean sections, and obstetric emergencies.
- Pain Management Clinics: organizes injectable medications and supplies for nerve blocks and interventional pain procedures.



Standards

The ToronCAB™ 1300 is manufactured under internationally recognized quality and medical-device management standards, and the manufacturer's certification scope covers medical trolleys.

- ISO 13485:2016 - Medical Devices: Quality Management Systems, Requirements for Regulatory Purposes

FEATURES

Anesthesia Cart - ToronCAB™ 1300 Key Features

01

Premium Casters

High-grade silent casters with 125mm diameter, smooth running and anti-winding.



02

Quiet Drawer Slide

Drawer with three-steps noise-less slide, smooth and quiet, with self-locking function.



03

Classified Storage

The trolley is equipped with 5 drawers as standard. The drawers have built-in 18 medicine units, and the size is adjustable for convenient medicine classification.



04

Side Working Table

The side pull-out working table extends the workspace, providing more convenience without increasing the volume of the trolley.



05

Structure

The interior of the car body is made of steel structure with firm welding and strong stability.



06

Safety Bumper

HDPE blow molding technology process, the material has strong toughness, so the whole trolley has excellent



The ToronCAB™ 1300 brings anesthesia-specific organization, durable construction, and on-board power together in a cart built for the operating room.

- Five Color-Coded Drawers: color-coded drawers assign drug classes to set positions, so the team reaches the right group by sight during a case.
- Sixteen Adjustable Compartments: built-in dividers split the drawers into resizable compartments, so facilities sort the formulary and adjust each space to fit.
- Upper Medicine Shelf: a back-mounted shelf holds ampoules, syringes, and airway supplies in clear bins, visible and within reach above the work surface.
- Side-Pull Work Table: a pull-out table adds preparation space and tucks away to keep the cart compact in the operating room.
- Central Control Lock: one lock secures all five drawers in a single action, keeping controlled and high-alert medications protected between cases.
- Built-In Power Strip: a side power strip runs monitors, pumps, or warmers at the cart without relying on wall outlets.
- Impact-Resistant HDPE Body: a blow-molded body resists knocks and wipes down easily, with a welded steel internal frame for rigidity.
- Silent Self-Locking Slides and 125 mm Casters: three-step self-locking drawer slides and high-grade silent casters keep operation smooth and quiet.
- Needle Disposal and Dual Waste Bins: an on-board sharps container and two waste bins handle sharps and waste safely at the point of use.

THEORY & METHOD

Theory and Method

Anesthesia care depends on selecting the right drug quickly and correctly, so the ToronCAB™ 1300 organizes medications by color and compartment. Five color-coded drawers assign drug classes to set positions, and each drawer divides into sixteen

adjustable compartments that a facility resizes to fit its formulary. Staff reach the correct group by color and find the exact drug in its compartment, rather than reading every label during induction. Organizing storage this way reduces hesitation and selection errors when timing matters.

Beyond the drawers, the cart adds organization and work space for a busy case. An upper medicine shelf carries ampoules, syringes, and airway supplies in clear bins, so the team sees and grabs small items at a glance above the work surface. A side-pull work table slides out to give extra preparation space, then tucks away so the cart stays compact in the operating room. Keeping supplies visible and work space ready helps the team move without interruption.

Build quality, security, and power round out the design. An HDPE blow-molded body gives the cart strong impact resistance and smooth, wipeable surfaces, while a welded steel internal frame holds everything rigid. A central control lock secures all drawers at once, so the team closes down controlled and high-alert medications whenever the cart is unattended. A side power strip runs monitors or pumps, and four 125 mm silent casters with brakes move the cart quietly and lock it firmly during a case.

TECHNICAL SPECIFICATIONS

Anesthesia Cart - ToronCAB™ 1300 Technical Specifications

Parameter	Specification
External Dimensions	L815 × W540 × H980 mm (±10 mm)
Body Construction	Blow-molded HDPE with welded steel internal frame
Medication Storage	Five color-coded drawers with sixteen adjustable compartments
Supply Organization	Upper medicine shelf with clear bins
Work Surface	Side-pull work table

Parameter	Specification
Locking	Central control lock for the drawers
Power	Side-mounted power strip
Waste and Sharps	Two waste bins plus a needle disposal container
Drawer Slides	Three-step, self-locking silent slides
Casters	Four × 125 mm silent swivel casters with brakes

Standard Configuration

Item	Quantity	Description
Drawer	5 pcs	Color-coded drawers with sixteen adjustable compartments
Central Lock	1 pc	Central control lock securing all drawers from one point
Sliding Working Table	1 pc	Side-pull work surface for drawing up drugs
Medicine Cassette Holder	1 set	Upper shelf with clear bins for ampoules, syringes, and supplies
Needle Disposal Container	1 pc	On-board sharps container for used needles
Waste Bin	2 pcs	Side-mounted bins for general and clinical waste
Power Strip	1 pc	Side-mounted strip powering monitors or devices at the cart
Individual Brake Caster	4 pcs	125 mm silent swivel casters, each with a brake

— • Complete Accessories • —

Thoughtful design for greater convenience and efficiency



 **Side Pull-out Working Table**

Hidden type working table, increasing the working area for nurse.



 **Power Strip**

Hidden type working table, increasing the working area for nurse.



 **Needle Disposal Container**

Safe disposal of sharps to help prevent cross-infection.



 **Waste Bins (2 pcs)**

Two separate waste bins for efficient waste classification and management.



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

