



TORONTECH



UTILITY CART (STAINLESS STEEL)

Utility Cart (Stainless Steel) - ToronCAB™ 1530

ToronCAB™ 1530 is a two-tier stainless steel utility cart for moving instruments, supplies, and equipment around wards, clinics, and sterile areas. Full stainless steel construction, including both decks, gives a hygienic surface that withstands repeated cleaning and disinfection. Guardrails around each deck keep items from sliding off, and a side push handle gives staff a firm grip when steering. Four 100 mm casters, two with brakes, roll the cart smoothly and hold it steady when parked.

Hospitals and clinics choose the Torontech ToronCAB™ 1530 where an all-metal cart suits the environment better than plastic. Stainless steel resists corrosion, tolerates strong disinfectants, and stands up to years of heavy use, which makes it a good fit for sterile and high-traffic areas. Corner bumpers on the frame cushion contact with walls and equipment during transport. The cart ships in knock-down form to lower freight cost, and open shelving lets it adapt to whatever a department needs to move.

The compact 730 by 470 mm footprint moves easily through corridors and doorways and fits into tight treatment, storage, or procedure areas. Built for daily use, the ToronCAB™ 1530 gives any department a durable, easy-to-clean cart for organized, mobile storage.

APPLICATIONS

Utility Cart (Stainless Steel) - ToronCAB™ 1530 Applications

The ToronCAB™ 1530 suits departments that need a durable, fully cleanable cart for moving instruments, supplies, or equipment. Its all-stainless build fits clinical and sterile settings where metal is preferred over plastic.

- Operating Rooms and Sterile Areas: transports sterile packs and instruments in an all-metal cart that tolerates strong disinfection.

- Central Sterile and Supply Departments: moves instrument sets through cleaning, packing, and distribution where hygiene standards are high.
- Treatment and Examination Rooms: carries instruments, dressings, and consumables to the bedside and between rooms.
- Wards and Nursing Stations: moves supplies, linens, or equipment during rounds and restocking.
- Laboratories: carries samples, glassware, and equipment between benches and storage.
- Dental and Outpatient Clinics: keeps tools and materials organized and mobile during procedures.
- Pharmacies and Storerooms: relocates stock and supplies during picking and restocking.



FEATURES

Utility Cart (Stainless Steel) - ToronCAB™ 1530 Key Features

The ToronCAB™ 1530 combines durable all-stainless construction, secure two-tier shelving, and easy mobility in a simple utility cart.

- Full Stainless Steel Construction: frame and both decks are stainless steel, giving a corrosion-resistant surface that tolerates strong disinfection.
- Two Decks with Guardrails: two open shelves with rails on each keep items organized and stop them sliding off in transit.
- Side Push Handle: a lateral handle gives staff a firm grip and control when steering the loaded cart.
- Corner Bumpers: bumpers on each corner cushion contact with walls and equipment during transport.
- Four 100 mm Casters with Brakes: four casters roll smoothly over hospital floors, and two with brakes hold the cart steady when parked.
- Knock-Down Construction: the cart ships flat to lower freight cost and assembles on site without special tools.

THEORY & METHOD

Theory and Method


The ToronCAB™ 1530 keeps transport simple by separating loads across two open decks rather than enclosing them in drawers. Two tiers double the usable surface within a small footprint, so staff group items by deck and load or unload from any side. A guardrail runs around each deck to stop items sliding off as the cart moves over thresholds and around corners. A side push handle gives the operator a fixed grip and good control when steering the loaded cart.

Material choice sets this cart apart. Full stainless steel construction, frame and decks alike, gives a single cleanable surface with no plastic to crack or stain, which suits sterile and frequently disinfected areas. Corner bumpers cushion knocks so the cart and the building both avoid damage, and four 100 mm casters carry it over hospital floors with brakes on two for stability. The frame assembles from a knock-down kit, which ships flat to cut freight cost and goes together on site.

TECHNICAL SPECIFICATIONS

Utility Cart (Stainless Steel) - ToronCAB™ 1530 Technical Specifications

Parameter	Specification
External Dimensions	730 (L) × 470 (W) × 975 (H) mm (±10 mm)
Construction	Full stainless steel (frame and decks)
Decks	Two stainless steel decks with guardrails
Handle	Lateral push handle
Corner Protection	Bumpers on each corner
Casters	Four × 100 mm; two with brakes
Construction Type	Knock-down (flat-pack)



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 

Phone: +971 488 19252
Email: me@torontech.com

