



TORONTECH

PULMONARY AEROSOLIZER



Pulmonary aerosolizers and dry powder insufflators are devices designed to deliver precise and measured aerosol samples directly into the trachea of rodents, enabling accurate drug delivery to the lungs. The operation table securely holds rats and mice in various positions to enhance procedure efficiency. The rodent laryngoscope allows for quick, minimally invasive insertion into the mouth and trachea, with controlled insertion depth. These tools are suitable for quantitative lung drug delivery in mice, rats, and guinea pigs.

APPLICATIONS

Pulmonary Aerosolizer
Pulmonary Aerosolizer

Pulmonary Aerosolizer
Pulmonary Dry Powder Insufflator

Pulmonary Aerosolizer

Pulmonary Aerosolizer

Application

- 1. Inhalation Therapy
- 2. Drug Research

3. Medical Protection

4. Biological and Chemical Hazard Assessment



Torontech Inc.

251 Consumers Road, Suite 1200
Toronto, Ontario, M2J 4R3, Canada 🇨🇦

Phone: +1 416 368 2721 | Fax: +1 416 981 7652
Email: info@torontech.com

Torontech, LLC

601 Brickell Key Drive, Suite 700
Miami, FL 33131, USA 🇺🇸

Toll-Free: 1-866-383-7919
Email: sales@torontech.com

Torontech ME FZE

Meydan Grandstand, 6th Floor,
Meydan Road, Nad Al Sheba,
Dubai-261440, UAE 🇦🇪

Phone: +971 488 19252
Email: me@torontech.com

